de

# Anti-A, Anti-B, Anti-AB

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CE 0483 IVD

# Anti-A, Anti-B, Anti-AB

Monoclonal test reagents for blood-group determination

INFORMATION FOR PROFESSIONAL USE

## Application

The monoclonal test reagents are used to detect the antigens of the ABO blood group system by agglutination of human red blood cells.

## Principle of the test

The human ABO blood group system is defined by presence or absence of the antigens A and/or B on the erythrocytes. The special feature of the ABO blood group system is that persons lacking the A and/or B antigens from the red blood cells regularly have antibodies in the serum directed against the missing antigen or antigens. The following table shows the principal antigens and antibodies of the ABO system:

Blood group	Antigens present on the red blood cells	Antibodies present in the serum
0	neither A nor B	Anti-A and Anti-B
A	A	Anti-B
B	B	Anti-A
AB	A and B	none

The specific monoclonal antibodies of the test reagents agglutinate red blood cells possessing the homologous antigen or antigens. No agglutination indicates the absence of the relevant antigen. The results of the identification of the antigens using the test reagents must be confirmed by testing the serum with red blood cells of a known ABO blood group (see „Interpretation of test results“ below). Any difference in the results of the red blood cell antigen and the alloagglutinin test must be clarified.

## Composition

The test reagents are manufactured from the cell culture supernatant of hybridoma cell lines. The hybridoma cell lines were obtained by fusion of spleen cells of immunized mice with mouse myeloma cells. The antibodies belong to the IgM isotype. The monoclonal test reagents contain antibodies of the following hybridoma cell lines: Test reagent Anti-A: clone sifin A-11H5, test reagent Anti-B: clone sifin B-6F9, test reagent Anti-AB: clone sifin A-5E10 and clone sifin B-2D7.

The test reagents Anti-A and Anti-B are coloured to avoid confusion and to permit checking of the test mixture through the blue and yellow colours, respectively.

Preservative: sodium azide ( $\leq 0.99$  mg/ml).

## Shelf life and storage

If stored at 2...8 °C, the test reagents can be used until the date indicated on the label. After opening, they should be stored well closed at 2...8 °C. While performing the test, the reagent can be stored for 2-3 h at 18...26 °C. Avoid bacterial contamination. Avoid bacterial contamination. Do not use if markedly turbid.

## Test material

Blood should be drawn by an aseptically, with or without an anticoagulant (EDTA). The specimen should be tested as soon as possible after collection. If testing is delayed, the specimen must be stored at 2...8 °C. Bacterial contamination of the specimen may cause false test results. Blood drawn into heparin should not be stored for longer than two days. Blood drawn into EDTA as well as clotted blood may be tested up to 14 days after collection, and donor blood may be tested up to the expiration date. Note, however, that storage may result in weaker than normal reactions (see „Notes“ below).

## Materials and equipment not supplied

Test tubes (75 × 12 mm or 70 × 8 mm), plates, pipettes (drop volume 40-50 µl), slides, applicator sticks, physiological sodium chloride solution (154 mmol/l), laboratory centrifuge, timer.

## Methods

### 1. Slide or plate test

- Place 1 drop (40-50 µl) of test reagents Anti-A, Anti-B and Anti-AB, respectively, on a properly labelled plate. Instead of the plate, slides (Lauer technique) may be used.
- To each drop of test reagent, add 1 drop (40-50 µl) of a 5% suspension of the red blood cells to be tested (in physiological saline). Mix carefully with a clean glass or plastic applicator stick.
- Read for agglutination after 5, at the most 20 min incubation at 18...26 °C by rotating or rocking the plate or slide slowly. Instead of the 5% suspension, a drop of blood may be used. Caution: Reaction times longer than 5 minutes may cause false positive results due to drying of the reactants.
- Record test results. Test results must be interpreted immediately upon completion of the test as the reaction product is unstable.

### 2. Tube test

- Prepare a 3-5% suspension of the red blood cells to be tested (washed or unwashed) in physiological sodium chloride solution.
- Place 1 drop (40-50 µl) of test reagents Anti-A, Anti-B and Anti-AB, respectively, in three small (e. g., 70 × 8 mm), properly labelled test tubes.
- Add 1 drop (40-50 µl) of the red blood cell suspension to each tube and mix thoroughly by shaking.
- Centrifuge immediately for 1 to 2 min at rcf 180 to 270 × g.
- Examine the supernatant for the absence of haemolysis. Haemolysis may be the consequence of bacterial contamination.
- Gently shake the tube to loosen the sediment from the bottom of the tube and examine for agglutination.
- Record test results. Test results must be interpreted immediately upon completion of the test as the reaction product is unstable.

## Interpretation of the test results

There must be agreement between the results of the antigen detection (cell grouping) and the detection of the alloantibodies (serum grouping) of a blood specimen. The following table shows the reaction patterns of the most common ABO phenotypes and their approximate frequencies:

Cell Grouping			Serum Grouping				ABO Group	Frequency %
Anti-A	Anti-B	Anti-AB	A <sub>1</sub> cells	A <sub>2</sub> cells	B cells	0 cells		
+	-	+	-	-	+	-	A	40
-	+	+	+	+	-	-	B	11
-	-	-	+	+	+	-	0	45
+	+	+	-	-	-	-	AB	4

The interpretation of reactions obtained when testing infant bloods may be complicated by the fact that the infant's serum does not necessarily contain antibody for any antigen absent from the cells, and passive anti-A and/or anti-B antibodies from the mother's circulation may yield conflicting results when tests are performed on umbilical cord blood specimens. Umbilical cord blood specimens may also produce weaker than normal reactions in the cell grouping, as the ABH antigens are imperfectly developed at birth.

## Warnings

- The use of biotechnological production means that the risk of contamination of the test reagents with infectious agents is virtually excluded. Since the test reagents contain materials of animal origin (fetal calf serum, stabilizer) they should be treated as a potential health hazard and be handled accordingly.
- The test reagents contain sodium azide; avoid any contact with the skin or mucous membranes.

## Notes and limits of the procedure

The performance assessment according to the Common Technical Specifications (CTS: decision of the EU commission on 03.02.2009) yielded the following calculated diagnostic values:

reagent	positive	false negative	sensitivity	negative	false positive	specificity
Anti-A	12359	0	100 %	13326	0	100 %
Anti-B	5535	0	100 %	20150	0	100 %
Anti-AB	1215	0	100 %	462	0	100 %

- The monoclonal test reagents agglutinate weakly expressed antigens with normal or weaker agglutination strength (A<sub>3</sub>, B<sup>weak</sup>) or weaker to negative reaction (A<sub>x</sub>). To determine weakly expressed antigens the tube test should be used because of its higher sensitivity, if necessary with incubation for 30 minutes.
- Monoclonal test reagent Anti-B does not react with acquired B.
- Monoclonal test reagent Anti-A reacts with Tn. Tn positive persons must be excluded from donating blood as the occurrence of Tn is considered to be a symptom of a preleukaemic state and the red blood cells are polyagglutinable.
- Titre determinations: To prepare dilutions of the monoclonal test reagents, the use of physiological saline containing 2 % protein (AB serum, bovine serum albumin) is recommended to facilitate the resuspension of the reaction products from the plate or the bottom of the tube.

## Notes

- Control:** For each test, positive and negative control red blood cells have to be tested in parallel.
- Concentrations of red blood cells that are lower than 3 % may cause weaker reactions on plates or slides.
- Older blood specimens may show weaker reactions than fresh specimens.
- On plates, blood samples may occasionally react by rouleaux formation, which can be mistaken for a weak agglutination and may incorrectly be read as a positive result. This phenomenon has non-immunological causes. Rouleaux formation also occurs in heparin blood, blood from patients treated with plasma expanders (e. g., dextran or hydroxyethyl starch) as well as patients with plasmacytoma (high protein content, changes in protein composition), oncological disease (pathological blood count) or coagulation dysfunctions. For these patients, use the tube test, since this phenomenon is usually not detectable in suspensions.
- The concentrations of red blood cell suspensions and the reaction times and conditions specified above must be observed in order to obtain correct results.

## Literature

(1) Richtlinie zur Gewinnung von Blut und Blutbestandteilen und zur Anwendung von Blutprodukten (Richtlinie Hämotherapie), current version (2) Technical Manual, 17th ed. Bethesda, American Assoc. of Blood Banks, 2011 (3) Messeter L et al: Mouse monoclonal antibodies with Anti-A, Anti-B and Anti-AB specificities, some superior to human polyclonal ABO reagents. Vox Sang 1984; 46, 185-194

## Package size

10 ml

## Explanation of the symbols used

Batch code (Lot)

Catalogue number

In Vitro Diagnostic Medical Device

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