NADAL® Mononucleosis Test (test cassette)

REF 252017N-05

 ¡Ordene aquí!

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15. Referenzen

1. Hickey SM, Strasburger VC. *What Every Pediatrician Should Know About Infectious Mononucleosis In Adolescents*. *Pediatr Clin North Am*. 1997;44(6):1541-56.
2. Omori M. Mononucleosis <http://www.emedicine.com/EMERG/topic309.htm>
3. Linde A. Diagnosis of Epstein-Barr-Virus-related diseases. *Scand J Infect Dis Suppl*. 1996; 100:83-8.
4. Papesch M, Watkins R. *Epstein-Barr-Virus infectious mononucleosis*. *Clin Otolaryngol*. 2001; 26(1):3-8.
5. CDC National Center for Infectious Diseases. EBV & IM: <http://www.cdc.gov/ncidod/diseases/ebv.htm>

Rev. 1, 2019-07-04 OM/UJa

1. Intended Use

The NADAL® Mononucleosis Test is a rapid, chromatographic immunoassay for the qualitative detection of heterophile antibodies to Epstein-Barr virus (EBV) in human whole blood, serum or plasma. The test is intended for use as an aid in the diagnosis of infectious mononucleosis (IM) and designed for professional use only.

2. Introduction and Clinical Significance

Infectious Mononucleosis (IM) is caused by the Epstein-Barr virus (EBV), which is a member of the herpesvirus family. The symptoms of IM are: fever, sore throat and swollen lymph glands. In very rare cases, heart or central nervous system problems may occur. Diagnosis of IM is made based on the presence of heterophile anti-EBV antibodies. Heterophile anti-EBV antibodies belong to the IgM class. They are present in 80-90% of acute IM cases and can be detected in 60-70% of patients during the first week of clinical illness.^{1,2,3,4}

The NADAL® Mononucleosis Test is a simple test that utilises an extract of bovine erythrocytes to qualitatively and selectively detect heterophile anti-EBV antibodies in whole blood, serum or plasma in just minutes.

3. Test Principle

The NADAL® Mononucleosis Test is a lateral flow immunoassay for the qualitative detection of heterophile antibodies to Epstein-Barr virus (EBV) in whole blood, serum or plasma.

In this test, antigens extracted from bovine erythrocytes are immobilised in the test line region of the test cassette. During testing, the specimen reacts with further antigens extracted from bovine erythrocytes and precoated onto particles which have been applied to the conjugate pad. The mixture migrates chromatographically along the membrane and interacts with the immobilised antigens. If the specimen contains heterophile anti-EBV antibodies, a coloured line will appear in the test line region, indicating a positive result. If the specimen contains no heterophile anti-EBV antibodies, no coloured line will appear in the test line region, indicating a negative result.

To serve as a procedural control, a coloured line should always appear in the control line region (C), indicating that a proper volume of specimen has been added and membrane wicking has occurred.

4. Reagents and Materials Supplied

- 5 NADAL® Mononucleosis test cassettes incl. disposable pipettes
- 1 buffer
- 1 positive control (diluted human plasma containing heterophile anti-EBV antibodies, 0.09% sodium azide)
- 1 negative control (diluted human plasma, 0.09% sodium azide)
- 1 package insert

5. Additional Materials Required

- Specimen collection containers
- Centrifuge (for serum/plasma)
- Timer
- Lancet (for fingerstick whole blood only)

- If necessary, heparinised capillary tubes and dispensing bulb (for fingerstick whole blood only)

6. Storage & Stability

Test kits should be stored at 2-30°C and used before the expiry date printed on the packaging. Test cassettes are stable until the expiration date printed on foil pouches. Test cassettes should remain in sealed pouches until use. Do not freeze tests. Do not use tests beyond the expiration date.

7. Warnings and Precautions

- For professional *in-vitro* diagnostic use only.
- Carefully read through the test procedure prior to testing.
- Do not use the test beyond the expiration date indicated on the package.
- Do not use the test if the foil pouch is damaged.
- Do not reuse tests.
- Do not add samples to the reaction area (result area).
- In order to avoid contamination, do not touch the reaction area (result area).
- Avoid cross-contamination of specimens by using a new specimen collection container for each specimen obtained.
- Do not eat, drink or smoke in the area where specimens and test kits are handled.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are being assayed.
- Handle all specimens as if they contain infectious agents. Observe established precautions for microbiological risks throughout all procedures and standard guidelines for the appropriate disposal of specimens.
- The test kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not completely guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious, and handled in accordance with usual safety precautions (e.g., do not ingest or inhale).
- Humidity and temperature can adversely affect test results.
- Used testing materials should be discarded according to local regulations.

8. Specimen Collection and Preparation

The NADAL® Mononucleosis Test can be performed using whole blood (from venipuncture or fingerstick), serum or plasma.

To collect venipuncture whole blood specimens:

Collect blood specimens with anticoagulants (sodium oxalate or lithium heparin, potassium or sodium EDTA, sodium oxalate, sodium citrate) following standard laboratory procedures.

To collect fingerstick whole blood specimens:

- Wash the patient's hand with soap and warm water or clean it with an alcohol swab. Allow it to dry.
- Massage the hand without touching the intended puncture site by rubbing down the hand towards the fingertip of the middle or ring finger.
- Puncture the skin with a sterile lancet. Wipe away the first sign of blood.
- Gently rub the hand from wrist to palm to finger to form a rounded drop of blood over the puncture site.

Separate serum or plasma from blood as soon as possible to avoid haemolysis. Use only clear, non-haemolysed specimens.

Testing should be performed immediately after specimen collection. Do not leave specimens at room temperature for prolonged periods of time. Serum and plasma specimens can be stored at 2-8°C for up to 3 days. For long-term storage, specimens should be kept below -20°C. Whole blood collected by venipuncture should be stored at 2-8°C if the test is to be run within 2 days of collection. Do not freeze whole blood specimens. Fingertick whole blood should be tested immediately. Bring specimens to room temperature prior to testing. Frozen specimens should be completely thawed and mixed well prior to testing. Specimens should not be frozen and thawed repeatedly. If specimens are to be shipped, they should be packed in compliance with federal regulations for transportation of etiologic agents.

9. Test Procedure

Bring tests, buffer, specimens and/or controls to room temperature (15-30°C) prior to testing.

1. Remove the test cassette from the sealed foil pouch and place it on a clean and level surface. The test cassette should be used immediately after opening the foil pouch. Label the test cassette with the patient or control identification.

2. a) For serum or plasma specimen:

Holding the dropper vertically, add 1 drop (approximately 25 µL) of serum or plasma to the sample well (S) of the test cassette.

b) For venipuncture whole blood specimen:

Holding the dropper vertically, add 2 drops (approximately 50 µL) of whole blood to the sample well (S) of the test cassette.

c) For fingerstick whole blood specimen:

Apply the end of the capillary tube to the blood until it is filled with approximately 50 µL. Avoid air bubbles. Place the bulb onto the top end of the capillary tube and squeeze the bulb to dispense the whole blood specimen into the sample well (S) of the test cassette.

3. Holding the buffer bottle vertically, add 1 drop (approximately 55 µL) of buffer to the sample well (S).
4. Start the timer.
5. Wait for the coloured line(s) to appear. Read the result after 5 minutes. Do not interpret the result after more than 10 minutes.

10. Result Interpretation

Positive:

Two coloured lines appear on the membrane. One line appears in the control line region (C) and the other line appears in the test line region (T).

Note:

The colour intensity in the test line region (T) may vary depending on the concentration of anti-EBV antibodies present in the specimen. Therefore, any shade of colour in the test line region (T) should be considered positive.

Negative:

One coloured line appears in the control line region (C). No apparent coloured line appears in the test line region (T).

Invalid:

The control line (C) fails to appear. Results from any test which has not produced a control line at the specified reading time must be discarded. Please review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit immediately and contact your distributor.

Insufficient specimen volume, incorrect operating procedure or expired tests are the most likely reasons for the control line failure.

11. Quality Control

An internal procedural control is included in the test cassette: A coloured line appearing in the control line region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

Good laboratory practice (GLP) recommends the use of control materials to ensure proper test kit performance.

Procedure for External Quality Control Testing

1. Holding each bottle vertically, add 1 drop (approximately 40 µL) of positive or negative control solution to the sample well (S) of the test cassette. Then, holding the bottle vertically, add 1 drop (approximately 55 µL) of buffer to the sample well (S).
2. Continue with step 4 of 'Test Procedure'.
3. If the controls do not yield the expected results, do not use the obtained results. Repeat the test or contact your distributor.

12. Limitations

- The NADAL® Mononucleosis Test is for *in-vitro* diagnostic use only. The test should only be used for the detection of anti-EBV antibodies in whole blood, serum or plasma specimens. Neither the quantitative value nor the rate of increase in concentration of anti-EBV antibodies can be determined by this qualitative test.
- The NADAL® Mononucleosis Test only detects the presence of anti-EBV antibodies in specimen and should not be used

as the sole criterion for diagnosis of infectious mononucleosis infection.

- As with all diagnostic tests, all results obtained with this test should be interpreted in conjunction with other clinical information available to the physician.
- If the test result is negative and clinical symptoms persist, additional testing using other clinical methods is recommended. A negative result does not at any time preclude the possibility of IM infection.

13. Expected Results

Infection with Epstein-Barr virus (EBV) during adolescence or young adulthood causes infectious mononucleosis (IM) in 35% to 50% of reported cases.^{1,5} The incidence of EBV-associated IM in the USA has been estimated at 45 per 100,000 and is the highest in adolescent and young adults - roughly 2 out of 1,000. No seasonal pattern of EBV infection exists. The incubation period is 10 to 60 days, though the incubation period of 7 to 14 days is common amongst children and adolescents.

14. Performance Characteristics

Sensitivity

The NADAL® Mononucleosis Test has been evaluated with specimens confirmed positive or negative by a leading, commercially-available slide agglutination test. The slide agglutination test served as a reference method for the NADAL® Mononucleosis Test. The result shows that the sensitivity of the NADAL® Mononucleosis Test is >99.9% relative to the slide agglutination test.

Specificity

The NADAL® Mononucleosis Test uses antigens that are highly specific for anti-EBV antibodies in whole blood, serum or plasma. The results show that the specificity of the NADAL® Mononucleosis Test is 98.9% relative to the slide agglutination test.

NADAL® Mononucleosis Test vs. slide agglutination test

NADAL® Mononucleosis Test	Slide agglutination test		
	Positive	Negative	Total
Positive	60	1	61
Negative	0	89	89
Total	60	90	150

Relative sensitivity: >99.9% (95.1% - 100%)*

Relative specificity: 98.9% (94% - 99.9%)*

Overall agreement: 99.3% (96.3% - 99.9%)*

*95% Confidence interval

Interfering substances

Negative and middle positive plasma and serum specimens (ELISA confirmed) were spiked with the analytes at the concentrations listed below and tested in triplicate.

Analyte	Concentration
Ascorbic acid	20 mg/ml
Haemoglobin	1000 mg/dl
Gentisic acid	20 mg/dl
Oxalic acid	60 mg/dl
Bilirubin	1000 mg/dl

Analyte	Concentration
Uric acid	20 mg/ml
Acetaminophen	20 mg/dl
Acetylsalicylic acid	20 mg/dl
Methanol	10%
Creatine	200 mg/dl
Albumin	2000 mg/dl
Caffeine	20 mg/dl

The analytes showed no interference with the NADAL® Mononucleosis Test. There were no obvious differences among the results obtained with three lots of tests.

Precision

Inter- and intra-lot variability

Three independent lots of the NADAL® Mononucleosis Tests were tested with negative as well as low, medium and high positive serum and plasma specimens in 10-fold determinations. All results matched the expectations. No deviating results were obtained within the 10-fold determinations or between the different lots.

Cross-reactivity

Samples, positive for the following parameters were tested with the NADAL® Mononucleosis Test:

RF	HCV
HbsAg	TB
HbeAg	HIV
HbeAb	Syphilis
HbcAb	

No cross-reactivity was observed indicating that the NADAL® Mononucleosis Test has a high degree of specificity for human anti-EBV antibodies.

15. References

1. Hickey SM, Strasburger VC. *What Every Pediatrician Should Know About Infectious Mononucleosis in Adolescents.* *Pediatr Clin North Am.* 1997;44(6):1541-56.
2. Omori M. Mononucleosis <http://www.emedicine.com/EMERG/topic309.htm>
3. Linde A. Diagnosis of Epstein-Barr-Virus-related diseases. *Scand J Infect Dis Suppl.* 1996; 100:83-8.
4. Papesch M, Watkins R. *Epstein-Barr-Virus infectious mononucleosis.* *Clin Otolaryngol.* 2001; 26(1):3-8.
5. CDC National Center for Infectious Diseases. EBV & IM: <http://www.cdc.gov/ncidod/diseases/ebv.htm>

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15. Lähteet

- Hickey SM, Strasburger VC. *What Every Pediatrician Should Know About Infectious Mononucleosis In Adolescents*. *Pediatr Clin North Am*. 1997;44(6):1541-56.
- Omori M. Mononucleosis <http://www.emedicine.com/EMERG/topic309.htm>
- Linde A. Diagnosis of Epstein-Barr-Virus-related diseases. *Scand J Infect Dis Suppl*. 1996; 100:83-8.
- Papesch M, Watkins R. *Epstein-Barr-Virus infectious mononucleosis*. *Clin Otolaryngol*. 2001; 26(1):3-8.
- CDC National Center for Infectious Diseases. EBV & IM: <http://www.cdc.gov/ncidod/diseases/ebv.htm>

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15. Referenser

1. Hickey SM, Strasburger VC. *What Every Pediatrician Should Know About Infectious Mononucleosis in Adolescents*. *Pediatr Clin North Am*. 1997;44(6):1541-56.
2. Omori M. Mononucleosis <http://www.emedicine.com/EMERG/topic309.htm>
3. Linde A. Diagnosis of Epstein-Barr-Virus-related diseases. *Scand J Infect Dis Suppl*. 1996; 100:83-8.
4. Papesch M, Watkins R. *Epstein-Barr-Virus infectious mononucleosis*. *Clin Otolaryngol*. 2001; 26(1):3-8.
5. CDC National Center for Infectious Diseases. EBV & IM: <http://www.cdc.gov/ncidod/diseases/ebv.htm>

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



15. Referanser

1. Hickey SM, Strasburger VC. *What Every Pediatrician Should Know About Infectious Mononucleosis In Adolescents.* *Pediatr Clin North Am.* 1997;44(6):1541-56.
2. Omori M. *Mononucleosis* <http://www.emedicine.com/EMERG/topic309.htm>
3. Linde A. *Diagnosis of Epstein-Barr-Virus-related diseases.* *Scand J Infect Dis Suppl.* 1996; 100:83-8.
4. Papesch M, Watkins R. *Epstein-Barr-Virus infectious mononucleosis.* *Clin Otolaryngol.* 2001; 26(1):3-8.
5. CDC National Center for Infectious Diseases. EBV & IM: <http://www.cdc.gov/ncidod/diseases/ebv.htm>

Rev. 1, 2019-07-04 KD

Symbol	Deutsch	English	Français	Español	Italiano	Polski
	CE Konformitätszeichen	CE marking of conformity	Conformité aux normes européennes	Conformidad europea	Conformità europea	Znak zgodności CE
	Gebrauchsanweisung beachten	Consult instructions for use	Consulter la notice d'utilisation	Consúltense las instrucciones de uso	Consultare le istruzioni per l'uso	Przestrzegać instrukcji obsługi
	In-vitro-Diagnostika	In-vitro diagnostic medical device	Dispositif médical de diagnostic in vitro	Producto sanitario para diagnóstico in vitro	Dispositivo medico-diagnostico in vitro	Tylko do diagnostyki in vitro
	Temperaturbegrenzung	Temperature limitation	Limites de température	Límite de temperatura	Limiti di temperatura	Temperatura przechowywania
	Chargenbezeichnung	Batch code	Code du lot	Código de lote	Codice lotto	Numer serii
	Nicht zur Wiederverwendung	Do not reuse	Ne pas réutiliser	No reutilizar	Non riutilizzare	Tylko do jednorazowego użytku
	Verwendbar bis	Use by	Utiliser jusqu'au	Fecha de caducidad	Utilizzare entro	Data ważności
	Bestellnummer	Catalogue Number	Référence du catalogue	Número de catálogo	Riferimento di Catalogo	Numer katalogowy
	Hersteller	Manufacturer	Fabricant	Fabricante	Fabbricante	Producent
	Ausreichend für <n> Ansätze	Sufficient for <n> tests	Suffisant pour "n" tests	Suficiente para <n> utilizaciones	Sufficiente per "n" saggi	Wystarczający na <n> Powtórzeń

Symbol	Português	Český	Suomi	Svenskt	Nederlands	Dansk	Norsk
	Conformidade com as normas europeias	CE certifikát	CE-merkitty	CE-märkning	CE-märkering	CE-mærkning	CE standardisert
	Consultar as instruções de utilização	Viz návod k použití	Katso käyttöohjetta	Läs bruksanvisningen	Raadpleeg de gebruiksaanwijzing	Se brugsanvisningen	Les bruksanvisning nøye
	Dispositivo médico para diagnóstico in vitro	Diagnostický zdravotnický prostředek in vitro	In vitro - diagnostikkaan tarkoitettu lääkinnällinen laite	Medicinteknisk produkt avsedd för in vitro-diagnostik	Medisch hulpmiddel voor in-vitrodiagnostiek	Medicinsk udstyr til in-vitro-diagnostik	In-vitro diagnostic medisinsk enhet
	Limites de temperatura	Teplotní omezení	Lämpötilarajat	Temperaturbegränsning	Temperatuurlimiet	Temperaturbegrænsning	Temperaturbegrensning
	Código do lote	Kód šarže	Eräkoodi	Satsnummer	Code van de partij	Batchkode	Merking
	Não reutilizar	Pro jednorázové použití	Kertäkäyttöinen	Får inte återvändas	Niet opnieuw gebruiken	Må ikke genbruges	Må ikke brukes om igjen
	Prazo de validade	Spotřebuje do	Käytettävä viimeistään	Används före	Houdbaar tot	Udløbsdato	Tidtaking
	Número de catálogo	Katalogov číslo	Luettelonumero	Listnummer	Catalogus nummer	Best il l ingsnummer	Katalog nummer
	Fabricante	Výrobce	Valmistaja	Tillverkare	Fabrikant	Fabrikant	Produsent
	Suficiente para <n> test	Dostačuje pro <n> testů	Lukumäärä <n> test	Räcker till <n> test	Voldoende voor <n> test	Tilstrækkeligt til <n> test	Tilstrekkelig for <n> tester

Our Teams

Germany:

Regensburg

Tel: +49 941 290 10-0
 Fax: +49 941 290 10-50

Moers

Tel: +49 2841 99820-0
 Fax: +49 2841 99820-1

Austria:

Tel: +49 941 290 10-29
 Free Tel: 0800 291 565
 Fax: +49 290 10-50
 Free Fax: 0800 298 197

UK & Ireland:

Tel: +49 941 290 10-18
 Free Tel – UK: 0808 234 1237
 Free Tel – IRE: 1800 555 080
 Fax: +49 290 10-50

France:

France Tel: 0800 915 240
 France Fax: 0800 909 493

Switzerland

Swiss Tel: 0800 564 720
 Swiss Fax: 0800 837 476

Belgium

Belgium Tel: 0800 718 82
 Belgium Fax: 0800 747 07

Luxembourg

Lux. Tel: 800 211 16
 Lux. Fax: 800 261 79

Spain:

Tel: +49 941 290 10-759
 Free Tel: 900 938 315
 Fax: +49 941 290 10-50
 Free Fax: 900 984 992

Italy:

Tel: +49 941 290 10-34
 Fax: +49 941 290 10-50

Poland:

Tel: +49 941 290 10-44
 Free Tel: 00 800 491 15 95
 Fax: +49 941 290 10-50
 Free Fax: 00 800 491 15 94

Portugal:

Tel: +49 941 290 10-735
 Tel. Verde: 800 849 230
 Fax: +49 941 290 10-50
 Fax Verde: 800 849 229

Netherlands:

Tel: +31 30 75 600
 Free Tel: 0800 0222 890
 Fax: +31 70 30 30 775
 Free Fax: 0800 024 9519

Nordic countries (Finland, Norway, Sweden, Denmark):

Tel: +31 703075 607
 Free Tel: +45 80 88 87 53
 Tax: +31 703030 775

Laboratory Diagnostics Team:

Tel: +49 941 290 10-40
 Fax: +49 941 290 10-50

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Carl-Zeiss-Strasse 12 • 47445 Moers • Germany

www.nal-vonminden.com • info@nal-vonminden.com

Fon: +49 2841 99820-0 • Fax: +49 2841 99820-1

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nal von minden GmbH