

QuikRead go® easy CRP



Remove the foil cover from the cuvette. Do not touch the clear flat surfaces in the lower part of the cuvette.



Pold QuikRead go Sample Collector 10μl almost horizontally just beneath the surface of the sample. Fill the capillary part of the sample collector completely.



Note! Check that the capillary part is completely filled. Do not wipe the sample collector! If there is excess sample on the sample collector, discard the sample collector.



Insert the sample collector into the cuvette within 2 minutes.



Close the cuvette tightly with the easy CRP reagent cap. Do not press down the turquoise inner part of the reagent cap. Perform the test within 2 hours.



Choose Measure on the display of the QuikRead go instrument. Insert the cuvette in to the measurement well of the instrument. The barcode should be facing you.



The result appears on the display and the cuvette rises up automatically when the measurement is completed.

Please consult the instructions for use before performing the test.

QuikRead go® easy CRP




Kit

Unopened kit at 2-25°C: until the expiry date of the kit



Cuvettes

Unopened foil pouch at 2-25°C: until the expiry date of the kit Unopened cuvettes after opening the foil pouch:





Opened cuvette

Let the buffer warm up to 18-25°C before testing Perform the test within 2 hours from opening the cuvette The sample is stable in the buffer for 2 hours



Reagent cap tube

Unopened tube at 2-25°C: until the expiry date of the kit

After the first opening at 2-25°C: 6 months



The sample collectors

Storage at 2-25°C: until the expiry date of

the Sample Collector 10µl box



CRP Controls

Unopened vial at 2-8°C: until the expiry date

Opened vial at 2-8°C: up to 2 months

Contact info

Aidian Oy

aidian@aidian.eu +358 10 309 3000

For orders: orders@aidian.eu

www.aidian.eu www.quikread.com

Product	Cat. no.
QuikRead go easy CRP with QuikRead go Sample Collector 10µl	153287
QuikRead go CRP Control Low	153765
QuikRead go CRP Control	153764
QuikRead go CRP Control High	153763
QuikRead go Instrument	133893
QuikRead go Sample Collector 10µl	153211

