# tamaVet<sup>®</sup>

## Parvovirus Test (test cassette)

REF. 1172033N-05 / REF. 1172633N-01

For veterinary use only

<b>DE</b>	Gebrauchsanweisung	2	<b>CZ</b>	Návod k použití	20
<b>EN</b>	Instruction for use	5	<b>NL</b>	Gebruiksaanwijzing	23
<b>FR</b>	Instructions d'utilisation	8	<b>PT</b>	Instruções de Utilização	26
<b>ES</b>	Instrucciones de uso	11	<b>TR</b>	Kullanma kılavuzu	29
<b>IT</b>	Istruzioni per l'uso	14		Symbols	35
<b>PL</b>	Sposób użycia	17		Our Teams	36

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Rev. 3, 2019-11-12 RZ

### 1. Intended Use

The tamaVet® Parvovirus Test is for use in detection of specific antigens of canine parvovirus and feline panleukopenia virus (CPV2a, 2b, 2c, as well CPV2 and FPV) **in dog and cat faeces**. The test is intended for use as an aid in the diagnosis of parvovirus (panleukopenia virus) infection.

### 2. Introduction and Clinical Significance

The canine parvovirus type 2 (CPV2) causes parvovirus, a highly contagious viral disease, which is particularly common amongst unvaccinated dogs. The virus can be fatal and can be found worldwide. The new types of CPV2 are also infectious for cats and can be transmitted between the two species. Feline panleukopeniavirus (FPV) is closely related to canine parvovirus and is a causative pathogen of feline panleukopenia (also known as feline distemper). This highly contagious disease is also found worldwide and mainly affects cats. CPV and FPV are very similar to each other and can be detected by the same rapid test as a consequence.

The clinical signs and the severity of the infection both in dogs and cats depend on the animal's age, breed and general health. The most severe infections occur in puppies and kittens.

The canine parvovirus particularly affect the gastrointestinal tract and the myocardium. Infected animals suffer from diarrhoea, vomiting and anorexia, as the gastrointestinal tract is affected. An infection with CPV2 can cause myocarditis in puppies younger than 8 weeks which can be fatal.

Infections with FPV are often asymptomatic and are particularly fatal in kittens between the age of 3 and 5 months. Infected animals suffer from diarrhoea, vomiting and anorexia, as it happens with CPV. This leads to severe dehydration.

The tamaVet® Parvovirus Test is a handy test cassette for the fast and on-site detection of CPV2a, 2b, 2c, as well CPV2 and FPV antigens.

### 3. Test Principle

The tamaVet® Parvovirus Test is a **sandwich immunoassay**. The assay employs specific antibodies which capture CPV/FPV-antigens from the sample and thus visualise the agent. The tamaVet® Parvovirus Test is a highly sensitive immunoassay and comprises a test strip housed in a handy test cassette.

The tamaVet® Parvovirus Test strip consists of a sample pad, a conjugate pad, a reaction membrane and an absorption pad. The conjugate pad contains specific gold-labeled anti-parvovirus antibodies. The reaction membrane is precoated with additional anti-CPV/FPV antibodies immobilised in the test line region and, as a control, antibodies immobilised in the control line region. The absorption pad takes up the fluid applied to the test cassette.

After application onto the test strip, CPV/FPV-antigens in the sample material react with the gold-labeled antibodies to form a complex in the conjugate pad. The fluid then migrates along the strip by capillary action. In the test line region, the formed complex is captured by the immobilised capture antibodies, causing the formation of a red line. The formation of a red line in the test line region indicates a positive result. If

the sample does not contain CPV/FPV-antigens, no line will form in the test line region. Additionally, a red line must form in the control line region irrespective of the presence of parvovirus antigens in the specimen. The control line serves as an internal control and indicates that sufficient sample fluid has been applied to the test strip and that membrane wicking has occurred.

### 4.1. Reagents and Materials Supplied (REF. 1172033N-05)

- 5 tamaVet® Parvovirus Test cassettes
- 5 sample tubes including an applicator stick with 1 mL buffer solution
- 5 cotton swabs
- 1 package insert

### 4.2. Reagents and Materials Supplied (REF. 1172633N-01)

- 1 tamaVet® Parvovirus Test cassette
- 1 sample tubes including an applicator stick with 1 mL buffer solution
- 1 cotton swabs
- 1 package insert

### 5. Additional Materials Required

Timer

### 6. Storage & Stability

The tamaVet® Parvovirus Test must be stored at room temperature or refrigerated (2-30°C). The test is stable until the expiration date printed on the foil pouch. The test cassette should remain in the sealed foil pouch until use. Do not freeze tests.

### 7. Warnings and Precautions

- For veterinary use only.
- For single use only.
- Store tests between 2-30°C. Do not freeze.
- Do not use tests beyond the expiry date.
- Do not use tests with damaged pouches.
- Perform the test within 60 minutes of removing the test cassette from the pouch.
- Consider the test results as invalid after more than 10 minutes.
- Please use the stated amount of sample. Too much faeces may interfere with the tests result. An incorrect number of drops can also lead to a false result.
- Place the test cassette on a smooth, flat surface while the test is being performed.
- Do not touch the reaction field.
- Use only the test components included in the kit. Do not mix or replace the test components with components from a different kit.
- Avoid cross contaminations by using a new sample tube and a new cotton swab for each sample.
- Sample materials are potentially infectious. During the test procedure, adhere to standard guidelines for the handling of potentially infectious material and chemical reagents. The use of protective clothing (laboratory coat, gloves, eye protection) is recommended. Materials that have come into

contact with the sample material should be disposed of in accordance with local regulations.

### 8. Specimen Collection and Preparation

The tamaVet® Parvovirus Test was developed for the examination of fresh **dog and cat faeces**.


#### General comments:

- Pooled faeces samples are suitable for testing.
- Mix the faecal sample thoroughly prior to testing.
- If the test is not to be performed on the day of sample collection, the faecal samples can be stored for a maximum of 2 days at 2°C to 8°C.
- The faecal sample must be brought to room temperature (15-30°C) before performing the test.
- Avoid picking up cat litter or other coarse particles during sampling.
- Use the cotton swab to take up faecal material. It is sufficient when the cotton tip is covered with faeces.

### 9. Test Procedure

1. Specimen and kit components must be brought to room temperature (15-30°C) before performing the test.
2. Open the aluminium pouch of the test cassette and place the cassette on a horizontal surface.
3. Open the sample tube and put the cotton swab with the faecal sample into the sample tube containing the reagent. Stir the fluid well with the cotton swab in order to dissolve the sample in the fluid. For an optimal test result, it is recommended to wring the swab into the sample tube by squeezing the tube with the swab in the middle several times to dissolve the sample. As a control, you should see the buffer solution discolor.
4. Remove the swab and close the sample tube. Shake the sample tube well. Break off the tip of the sample tube by pressing it firmly.
5. Holding the tube vertically with the lid down over the sample well, dispense **3 drops of the sample solution** to the sample well of the test cassette. Start the timer.

If the liquid does not migrate sufficiently along the test strip after a few seconds, add another drop to the sample well until the liquid starts running. If too much faeces is in the sample well you can press with the

upper part of the cotton swab into the sample well to reactivate the flow.

6. Read the result 10 minutes after starting the timer. **Do not interpret the result after more than 10 minutes.** The test and the result are considered as invalid after this period.

### 10. Result Interpretation

#### Positive result

The test is positive when the control line (C) and the test line (T) are visible in the reaction area. This indicates that CPV/FPV-antigens have been detected in the sample. If a weakly-defined test line appears, the test result is nevertheless positive.

#### Negative result

The test is negative when only the control line (C) is visible in the reaction area. No CPV/FPV-antigens have been detected in the sample material.

#### Invalid result

If no control line is visible after the test was performed, the test is invalid. All results that show no control line after the denoted time must be rejected. Review the test procedure and repeat the test with a new test cassette.

If problems continue, discontinue using the test kit and contact your distributor.

### 11. Quality Control

Several internal controls are included in the test cassette:

A red line appearing in the control line region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique. In addition, there is a green/blue line in the control line region visible before starting the test. It is considered a quality control method and will be washed away by the sample liquid during the test.

### 12. Limitations

- Test results must always be interpreted in the light of all the medical findings and never individually. All test results which do not correlate with the clinical picture must be repeated. Positive results detected with the tamaVet® veterinary rapid tests should be confirmed with another analytical method.
- Excretion of CPV/FPV is typically starts from day 4 to 7 after infection, which usually correlates with the onset of clinical symptoms. Virus excretion decreases from day 8 after infection. Therefore, it is important to collect faeces for viral detection at the beginning of clinical illness and if the test result is negative for parvovirus/panleukopenia virus, repeat the test after 1–2 days.

- There is a possibility that the test result can be distorted by technical errors, errors in the experiment or by substances or factors that influence the test and have not been mentioned here.

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nal von minden GmbH