

STAINLESS STEEL TRAY



DESCRIPTION

The **Teqler stainless steel tray** with higher edges is made of stainless steel 304 and can be both disinfected and sterilised. Available in different sizes.

PRODUCT DETAILS

- Material: stainless steel 304
- Rust-proof, even after regular reprocessing
- Higher edges to prevent instruments from rolling down
- Suitable for holding surgical instruments and medical consumables
- Can be disinfected and sterilised
- 1 pc.

CHOOSE YOUR SIZE

- 22.5 x 16 x 1.4 cm (LxWxH)
- 28.5 x 19 x 1.8 cm (LxWxH)
- 34 x 24.5 x 1.9 cm (LxWxH)
- 39 x 27.5 x 1.9 cm (LxWxH)

Overview




Article Number:
T134800-T134803

Category:
Stainless Steel Holloware

Dimensions:
Material: Edelstahl 304

Teqler is a trademark of NetMed S.à.r.l.

 NetMed S.à.r.l.
8, Am Scheerleck, 6868 Wecker, Luxembourg

 +352 267149 09  +352 267149 19  info@netmed.lu

Social Media

