

BEAKER GLASS, LOW FORM



DESCRIPTION

Stable, low-form glass beaker with high tolerance for sudden temperature changes, high temperatures and chemicals. With easily readable and detailed scale. Available in various sizes.

PRODUCT DETAILS

- Material: 3.3 borosilicate glass
- Detailed scale (white, non-smudging)
- High tolerance for heat and chemicals
- Low thermal expansion
- Crimped rim
- With spout
- Very stable
- Available in the following sizes: 100 ml, 250 ml, 1,000 ml




Overview

Article Number:
T135830-T135832

Category:
Lab Equipment

Teqler is a trademark of NetMed S.à.r.l.

 NetMed S.à.r.l.
8, Am Scheerleck, 6868 Wecker, Luxembourg

 +352 267149 09  +352 267149 19  info@netmed.lu

Social Media

