

SPECIFICATION

Labocult Laboratory Incubator



Compact laboratory incubator for practice and laboratory

The current regulation for laboratory incubators states that culture mediums may no longer stay in contact with the housing of the incubator. The Labocult has got from the beginning a shelf, which prevents the direct contact between culture mediums and housing.

Technical specifications

Voltage: 200 - 240 Volt

Frequency: 50 Hz

Power consumption: 300 watts max. / typical 80 Watt

Temperature code: max. 70 °C

Accuracy: $\pm 0,5$ °C at 37 °C

Timer: 0 - 9999 minutes

Internal dimensions: 230 x 200 x 200 mm

External dimensions: 300 x 330 x 330 mm

Weight: 9 kg

The Labocult heating oven features a digital control with timer shut-off function and easy-to-read LED display for simple temperature control. The display and control foil are washable, easy to use and include signal lamps for function and alarm in addition to the on/off switch. A new feature is the door with integrated viewing window and a button that can be reached from the outside to activate the interior lighting. This helpful function allows the inspection of the culture media without interrupting the incubation process. Temperatures, from 5 °C above room temperature to 70 °C, can be set in precise 0.1 °C increments. The housing of the Labocult is made of painted sheet steel, the interior is made of stainless chrome steel and has a usable depth of 20 cm. The digital control guarantees a short heating time and has an accuracy of ± 0.5 °C at a temperature of 37 °C.

Automatic temperature control!

The current regulation for laboratory incubators is stating that the breeding grounds may no longer be in direct contact with the bottom surface. The shelf prevents this contact.

WEEE Registrierungs-Nr. des Herstellers: DE74434443

scope of delivery

1 shelf for plates and immersion nutrient tray, power cord, operating instructions

SPECIFICATION

Labocult Laboratory Incubator

Type	PZN	REF	Unit
Labocult Laboratory incubator		H7 70	1 piece

