

# NEUATION



AMAZING  
PERFORMANCE IN  
A COMPACT PACKAGE

**iShak BL UNO VT**  
(Compact Microplate Shaker)

## iShak BL UNO VT (Compact Microplate Shaker)

Every microplate shaker advertises high performance and a long life, but look a little closer and you'll discover that these devices are far from what they promise. Low-quality components and low quality materials make these devices fragile and unreliable which is definitely not good news for the user. The cheap and low-quality motors used in these devices struggle with efficiency and power when they are stretched to their limits and more often than not they require constant maintenance to keep them running. Moreover, the low orbital diameter and the weak motor combine to provide satisfactory mixing of solutions leading to more problems for the user.


Neuation's expertise in understanding customers' expectations and requirements has led to the development of a wide-range of quality devices. iShak BL Uno VT is also one such device which has been designed to offer high-performance and exceptional ease-of-use in a compact package. The powerful brushless DC motor and variable parameter settings of iShak BL Uno VT make it an amazingly versatile device to suit all types of requirements.



## Product Features



Brushless DC Motor For Quiet Performance And A Maintenance-free Long Life.



Microprocessor Controlled



Variable Speed Setting Allows The User To Set Speeds From 300 To 1800 RPM



Variable Timer Setting Allows The User To Set The Timer From 1 To 99 Minutes Along With An Infinite Mode.

Orbital Diameter Of 2 MM



Digital Display



Inbuilt Counter Balance Ensures Stable Operation For Complete Peace-of-mind.



Allows Easy And Spill-free Loading/removal Of Plate With Turn Down Toggle. Clamp-tool Free Clamp Change.



## Product Features



Compact Footprint Saves Valuable Bench Space And Can Also Be Used Inside Cold Rooms/incubators (Upto 40° C)



Comes With Two Sets Of Interchangeable Clamps For Different Microplates

Safe Low Voltage Device



## Specifications

Specifications	iShak BL Uno VT
Speed Setting	300 to 1800 RPM
Orbital Diameter	2 mm
Motor type	BLDC
Timer setting	1 to 99 minutes
Digital Display	Yes
Inbuilt Counter balance	Yes
Dimension (W x D x H)	127 x 86 x 85 mm
Weight	1.36 kg
Attachments (included in standard park)	Replaceable clamps to use different types of standard micro plates
Included Attachment	Attachment for 20 x 1.5 / 2.0 ml micro tubes

## Ordering Information

iShak BL Uno VT		
Product Name	Cat #	Details
iShak BI UNO VT	iShak-BL-UNO-VT (Interchangeable Pin)	Single Microplate Shaker

Accessories		
Part Name	Cat #	Details
Attachment for Micro tubes	PS01-A1	20 slot Microtube Attachment
Clamps	PS01-C1	Clamps set for outside hold
Clamps	PS01-C2	Clamps set for inside hold

# NEUATION

[www.neuation.com](http://www.neuation.com)

AN ISO 9001 : 2015 & ISO 13485 : 2016 CERTIFIED MFG. COMPANY

\*We also execute OEM projects. Specifications can be changed without notice for quality improvement.