



**SUZHOU SUNING UNDERPAD CO., LTD**  
**蘇州市蘇寧床墊有限公司**

## Technical data sheet

<b>Item No.:</b>	<b>190649</b>
<b>Product name:</b>	Underpad
<b>Product description:</b>	Underpad for incontinence usage, mostly in hospital, senior homes, homecare and so on.
<b>Sizes of pad:</b>	60 x 90 cm
<b>Inner absorbent pad size:</b>	50 x 75 cm
<b>Raw Materials:</b>	Nonwoven, Tissue, PE, Fluff pulp, SAP
<b>Total weight per piece:</b>	90 g (28.7 g Amount of top nonwoven+tissue+ PE) (56.3 g Amount of fluff pulp) (5 g Amount of SAP) (22 gms Back layer PE) (13 gms Tissue) (12 gms Top layer nonwoven)
<b>Absorption:</b>	1.3 liters absorption capacity
<b>Color:</b>	Green
<b>Packaging:</b>	10 pcs / bag 10 bags / carton
<b>Packing Materials:</b>	<b>Polyethylene</b> <b>40 µm thickness</b>
<b>Carton size:</b>	42 x 29.5 x 46 cm
<b>Useful Life:</b>	<b>5 years</b>
<b>Handling and Storage:</b>	<b>Store dry, Avoid direct sunlight</b>
<b>Storage temperature:</b>	<b>temperature between -5-40 °C</b>

### Cautions, Warnings or contraindications

No. 70 Yong'an Road, Xuguan Industrial Park, High-Tech District, Suzhou, Jiangsu, China 215151  
Telephone: 86-512-66653331, Fax: 86-512-66651066  
[www.sn-anpa.com](http://www.sn-anpa.com)



**SUZHOU SUNING UNDERPAD CO., LTD**  
**蘇州市蘇寧床墊有限公司**

Not to use if the product is damaged or dirty, or do have other functional failures.

### **Rubbish and Final Disposition**

The product must be disposed, disposable waste regulation, after using.

**Date:** 2021-11-12

