


SURGICRYL[®]

MONOFAST

POLYGLECAPRONE 25



Type	Monofilament
Composition	Polyglecaprone 25, a copolymer made of 75% glycolide and 25% ε-caprolactone
Coating	None
Colour	Violet / Beige - undyed
Tissue reaction	Minimal
Absorption	The hydrolytic action by which the material is broken down results in total absorption in approximately 90 to 110 days. Approximately 75% of tensile strength remain after 7 days.
Presentation	Box with 12 sutures; Cassette

PRAXISDIENST
Medical Supplies since 1953
[Order here!](#)



*For quality
&
safety*

smi

www.sutures.be

SURGICRYL[®]

MONOFAST

POLYGLECAPRONE 25

CHARACTERISTICS

- High tensile strength
- Synthetic
- Breaks down by hydrolysis
- Monofilament with microscopic uniform diameter

- Pliability and softness

- Hermetically sealed packing

BENEFITS

- Excellent tissue support
- Minimal tissue reaction
- Predictable absorption rate
- Soft passage through the tissues without sawing, tissue drag and trauma
- Absence of capillarity (no support of bacterial growth)
- Excellent knot tie down

- Good handling properties
- Easy and secure knotting

- Guaranteed seal and product sterility


INDICATIONS Superficial soft tissue and/or ligation
 especially in General surgery
 Skin closure
 Plastic surgery
 Gastrointestinal surgery
 Gynaecology
 Urology

CONTRAINDICATIONS Cardiovascular surgery
 Neurosurgery
 Ophthalmic surgery


Being absorbable not to use when an extended wound support is required.


Caution in patients with delayed wound healing (e.g. compromised immune system).

TENSILE STRENGTH



ABSORPTION



 Without needle	Metric	USP	Length		Order number
	3	2/0	150 cm	<i>violet</i>	16300015
	3.5	0	150 cm	<i>violet</i>	16350015
	4	1	150 cm	<i>violet</i>	16400015
	5	2	120 cm	<i>violet</i>	16500012
	5	2	150 cm	<i>violet</i>	16500015

 Without needle	Metric	USP	Length		Order number
--	--------	-----	--------	--	--------------

Ligapack














5	2	250 cm	<i>violet</i>	16500025
---	---	--------	---------------	----------










 Round bodied taper point	Metric	USP	Length		Order number
--	--------	-----	--------	--	--------------









1/2 circle










HR 13 mm		0.7	6/0	45 cm	<i>violet</i>	16070113
		1	5/0	45 cm	<i>violet</i>	16100113
		1.5	4/0	75 cm	<i>violet</i>	16150113
HR 17 mm		1	5/0	45 cm	<i>violet</i>	16100117
		1.5	4/0	75 cm	<i>violet</i>	16150117
		2	3/0	75 cm	<i>violet</i>	16200117
HR 22 mm		1.5	4/0	75 cm	<i>violet</i>	16150122
		2	3/0	75 cm	<i>violet</i>	16200122
		3	2/0	75 cm	<i>violet</i>	16300122
HR 26 mm		1.5	4/0	75 cm	<i>violet</i>	16150126
		2	3/0	75 cm	<i>violet</i>	16200126
		3	2/0	75 cm	<i>violet</i>	16300126
		3.5	0	75 cm	<i>violet</i>	16350126
HR 30 mm		4	1	75 cm	<i>violet</i>	16400126
		2	3/0	75 cm	<i>violet</i>	16200130
		3	2/0	75 cm	<i>violet</i>	16300130
		3.5	0	75 cm	<i>violet</i>	16350130
HR 36 mm		4	1	75 cm	<i>violet</i>	16400130
		3	2/0	75 cm	<i>violet</i>	16300136
		3.5	0	90 cm	<i>violet</i>	16350136
HR 40 mm		4	1	90 cm	<i>violet</i>	16400136
		3	2/0	75 cm	<i>violet</i>	16300140
		3.5	0	90 cm	<i>violet</i>	16350140
	4	1	90 cm	<i>violet</i>	16400140	




 Round bodied taper point	Metric	USP	Length	Order number	
1/2 circle 	HR 48 mm 	3.5	0	90 cm	<i>violet</i> 16350148
		4	1	90 cm	<i>violet</i> 16400148
	HR 65 mm 	3	2/0	90 cm	<i>violet</i> 16300165





 Reverse cutting	Metric	USP	Length	Order number	
3/8 circle 	DS 12 mm 	0.7	6/0	45 cm	<i>undyed</i> 16071512
		0.7	6/0	75 cm	<i>undyed</i> 16071512-75
		0.7	6/0	75 cm	<i>violet</i> 16071512V
		1	5/0	45 cm	<i>undyed</i> 16101512
		1	5/0	45 cm	<i>violet</i> 16101512V
		DS 12 mm Easy Pass 	0.7	6/0	50 cm
		1	5/0	50 cm	<i>undyed</i> 16101512P
		1.5	4/0	50 cm	<i>undyed</i> 16151512P
	DS 16 mm 	1	5/0	75 cm	<i>undyed</i> 16101516
		1	5/0	75 cm	<i>violet</i> 16101516V
		1.5	4/0	75 cm	<i>undyed</i> 16151516
		1.5	4/0	75 cm	<i>violet</i> 16151516V
		2	3/0	75 cm	<i>undyed</i> 16201516
		DS 16 mm Easy Pass 	1	5/0	50 cm
		1.5	4/0	50 cm	<i>undyed</i> 16151516P
	DS 16 mm black 	1	5/0	75 cm	<i>violet</i> 16101516VB
	DS 19 mm 	1	5/0	45 cm	<i>undyed</i> 16101519
		1	5/0	75 cm	<i>violet</i> 16101519V
1.5		4/0	75 cm	<i>undyed</i> 16151519	
1.5		4/0	75 cm	<i>violet</i> 16151519V	
2		3/0	75 cm	<i>undyed</i> 16201519	
2		3/0	75 cm	<i>violet</i> 16201519V	
3		2/0	75 cm	<i>undyed</i> 16301519	
DS 19 mm Easy Pass 	1.5	4/0	50 cm	<i>undyed</i> 16151519P	
	2	3/0	50 cm	<i>undyed</i> 16201519P	

 Reverse cutting		Metric	USP	Length		Order number
	DS 19 mm black 	1	5/0	75 cm	<i>violet</i>	16101519VB
		1.5	4/0	75 cm	<i>violet</i>	16151519VB
	DS 24 mm 	2	3/0	75 cm	<i>undyed</i>	16201524
		3	2/0	75 cm	<i>undyed</i>	16301524
	DS 24 mm Easy Pass 	2	3/0	75 cm	<i>undyed</i>	16201524P-75
	DS 26 mm 	1.5	4/0	75 cm	<i>undyed</i>	16151526
		2	3/0	75 cm	<i>undyed</i>	16201526
		2	3/0	75 cm	<i>violet</i>	16201526V
		3	2/0	75 cm	<i>undyed</i>	16301526
	DS 30 mm 	3	2/0	75 cm	<i>violet</i>	16301526V
		3.5	0	75 cm	<i>violet</i>	16351530V
	DS 40 mm 	3.5	0	75 cm	<i>undyed</i>	16351540
4		1	90 cm	<i>undyed</i>	16401540	

 Reverse cutting		Metric	USP	Length		Order number
	HS 36 mm 	3	2/0	75 cm	<i>undyed</i>	16301436
		3.5	0	75 cm	<i>undyed</i>	16351436
	HS 40 mm 	4	1	90 cm	<i>undyed</i>	16401440

 Cutting		Metric	USP	Length		Order number
	GS 60 mm 	2	3/0	75 cm	<i>undyed</i>	16202160
		3	2/0	75 cm	<i>undyed</i>	16302160

 Round bodied taper cutting point		Metric	USP	Length		Order number
	DRT 16 mm 	1	5/0	75 cm	<i>violet</i>	16100416

 Blunt		Metric	USP	Length		Order number
1/2 circle 	HRB 40 mm		3.5	0	90 cm	<i>violet</i> 16351840
	HRB 48 mm		4	1	90 cm	<i>violet</i> 16401848

IN CASSETTES

Use within 3 months after first opening

Metric	USP	Length	Order number
1.5	4/0	15 m	16151
2	3/0	15 m	16201
2	3/0	25 m	16202
3	2/0	15 m	16301
3	2/0	25 m	16302
3.5	0	15 m	16351
3.5	0	25 m	16352
4	1	15 m	16401
4	1	25 m	16402
4	1	50 m	16405
5	2	15 m	16501
5	2	25 m	16502
5	2	50 m	16505
5	2	100 m	165010



PRAXISDIENST
Medical Supplies since 1953

Order here!