



# *TT-32DX Veterinary Dental X-ray Unit*

## **Characteristics and Usage**



This equipment is a high frequency dental X-ray machine designed and manufactured specially for pets, which has an attractive and tiny shape . Further more , it's pretty light and easy to carry for its usage.

The device works with the international DC and high-frequency power supply, all its components are installed on the central processing circuit board .

The product's used mainly to photograph the pet's intraoral tissues , the root canal treatment, among others . On the other hand, this machine can also be used in conjunction with sensors and image scanning equipment to greatly assist the users in the diagnosis and treatment process.



## TT-32DX Veterinary Dental X-ray Unit



Tube Voltage	60KV
Tube Current	2.0mA
Exposure Time	0.1-1.2s
Frequency	30KHz
Rated Power	120VA
Focus distance to skin	130mm
Tube Focus	0.4mm*0.4mm
Battery	2200mAh
Input Voltage	AC100V-240V±10%
Output Voltage	DC14.8V-16.8V
Ambient temperature	5°C-40°C
Relative Humidity	<75%RH
Net weight	2.2KG
Gross weight	4.5KG
Size(mm)	170*153*138
Package size(mm)	365*330*275