


Ambu® King Mask

Disposable anaesthesia face mask



Key Benefits

Market Leading Design

- Cushion is pliable and tacky to the touch, for an excellent face seal
- Ribbed on top of the cushion, for a no-slip grip
- Anatomically shaped cuff enabling a tight seal with minimum applied pressure
- Flexible dome
- Clear dome for easy observation of the patient's condition
- Latex free
- MR Safe
- Phthalate-free material
- Color-coded hook rings for easy size identification*

Ambu® King Mask

Because of the innovative quality design of the Ambu® King Mask, it is the market leading anaesthesia mask design available today. The Ambu® King Mask has a thin and tacky cushion which improves the ability to create a ventilation seal on a variety of patients. The design with rounded edges secures a smooth grip and user comfort.

The Ambu® King Mask is made of phthalate-free PVC material for patient and user safety. Furthermore, the Ambu® King Mask can be used in a MR environment, as it is MR Safe.

The Ambu® King Mask is also available in a scented version with Fresh Scent, Bubble Gum, Strawberry and Cherry scents. The Sweet Dreams™ masks release soothing aroma from natural oils in mask cushion. This can have a calming effect during induction, especially for children¹.

Five Scents

- Unscented
- Fresh Scent (mint)
- Bubble Gum
- Strawberry
- Cherry

Seven Sizes








- Size 1 Neonate
- Size 2 Infant
- Size 3 Toddler
- Size 4 Child / Small Adult
- Size 5 Medium Adult
- Size 6 Large Adult
- Size 7RD XL Adult

1. Carson KD et al. Quantitative assessment of patient preference for anaesthesia face mask. Irish Journal of Medical Science 1995;164 suppl. 14

* Not available for sizes 1 and 2.

Ambu[®] King Mask

Disposable anaesthesia face mask

Mask	Size 1 Neonate	Size 2 Infant	Size 3 Toddler	Size 4 Child / Small Adult	Size 5 Medium Adult	Size 6 Large Adult	Size 7RD XL Adult
Packaging Color Code							
Original (unscented)	1012	1022	1032	1045	1055	1065	1075
Fresh Scent (mint)	2012	2022	2032	2045	2055	2065	2075
Bubble Gum		K1122	1132	1145			
Strawberry		1222	1232	1245			
Cherry		1322	1332	1345			
Scented masks variety pack		1423	1433	1443			

Part #'s ending in 2 = 20 per case • Part #'s ending in 5 = 50 per case

Part #'s ending in 3 = 30 per case (Scented mask variety pack 10 of each of Bubble Gum, Strawberry, & Cherry)



Unscented size 5



Fresh Scent



Strawberry



Cherry



Bubble Gum

Materials

Part	Material
Cushion	PVC Polyvinyl Chloride (Phthalate-free)
Crown	Polycarbonate
Hook ring	Polycarbonate
Valve Body	Polycarbonate
Valve Stem	Silicone Rubber
Valve Plug	PP Polypropylene

Product weight

Size	Weight
1	15g
2	18.4g
3	32.4g
4	38.3g
5	50.6g
6	57.9g
7RD	51.3g

Ambu A/S

Baltorpbakken 13
DK - 2750 Ballerup
Denmark
Tel. +45 7225 2000
Fax +45 7225 2053
www.ambu.com



2797 US: Rx only